


## PE/Cyanine7 Rabbit Anti-Human FcεR1A Monoclonal Antibody For Fc - P12319

Item Number: CM0010229



### Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human FcεR1α (Fc-epsilon RI-alpha), the high-affinity IgE receptor. Validated for flow cytometry (FC). Ideal for allergy and mast cell/basophil research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE/Cyanine7 Rabbit anti-Human FcεR1α mAb                                                                            |
| Remarks/Alias    | FCE1A; FcERI; FCERIA                                                                                                |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 2205                                                                                                                |
| GeneID           | 2205                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-205 of human FcεR1α (NP_001992.1). |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | FC                                                                                                                  |
| Cross Reactivity | Human                                                                                                               |
| SWISS            | P12319                                                                                                              |
| Protein Weight   | 30kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |