

PE/Cyanine7 Rabbit Anti-Human Ifn-Gamma Monoclonal Antibody For Intracellular Flow Cytometry - P01579

Item Number: CM0005730



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human IFN-gamma. Validated for intracellular flow cytometry, ideal for detecting cytokine production in activated T and NK cells during immune response studies.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human IFN-gamma mAb
Remarks/Alias	IFG; IFI; IMD69
Species	Human
Gene ID (Human)	3458
Gene ID	3458
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN-gamma (NP_000610.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (Intracellular)
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag