

## Pe/Cyanine7 Rabbit Anti-Human Lilrb4/Cd85K/Ilt3 Mab - Q8Nhj6

Item Number: CM0012662



### Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human LILRB4 (CD85k/ILT3), a key inhibitory immune receptor. Validated for flow cytometry, suitable for studies on immune tolerance, monocyte regulation, and cancer immunology.

[Learn More](#)

| Parameter        | Value                                                                                                                            |
|------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE/Cyanine7 Rabbit anti-Human LILRB4/CD85k/ILT3 mAb                                                                              |
| Remarks/Alias    | B4; ILT3; LIR5; CD85K; ILT-3; LIR-5                                                                                              |
| Species          | Human                                                                                                                            |
| GeneID (Human)   | 11006                                                                                                                            |
| GeneID           | 11006                                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-259 of human LILRB4/CD85k/ILT3 (NP_001265355.2) |
| Source           | Rabbit                                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                                            |
| Application      | FC                                                                                                                               |
| Cross Reactivity | Human                                                                                                                            |
| SWISS            | Q8NHJ6                                                                                                                           |
| Protein Weight   | 49kDa                                                                                                                            |
| Shipping         | Ice bag                                                                                                                          |