

PE/Cyanine7 Rabbit Anti-Human Pd-L1/Cd274 Mab For Flow Cytometry - Q9Nzq7

Item Number: CM0009410



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human PD-L1/CD274 (Q9NZQ7). Designed for flow cytometry, this antibody aids in studying immune checkpoint pathways and tumor immune evasion in research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human PD-L1/CD274 mAb
Remarks/Alias	B7-H; B7H1; PDL1; PD-L1; hPD-L1; PDCD1L1; PDCD1LG1
Species	Human
GeneID	29126
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human PD-L1/CD274 (Q9NZQ7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NZQ7
Protein Weight	33kDa
Shipping	Ice bag