


Pe/Cyanine7 Rabbit Anti-Human Tim-3/Havcr2 Monoclonal Antibody For Flow Cytometry - Q8Tdq0

Item Number: CM0008561



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody specific for human TIM-3/HAVCR2 (CD366). This reagent targets the extracellular domain and is validated for flow cytometry, making it ideal for immune checkpoint analysis in cancer and chronic infections.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human TIM-3/HAVCR2 mAb
Remarks/Alias	TIM3; CD366; KIM-3; SPTCL; TIMD3; Tim-3; TIMD-3; HAVcr-2
Species	Human
Gene ID (Human)	84868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-200 of human TIM-3/HAVCR2 (NP_116171.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TDQ0
Protein Weight	33kDa
Shipping	Ice bag