


Pe/Cyanine7 Rabbit Anti-Human Tnf-A Monoclonal Antibody For Fc (Intra) - P01375

Item Number: CM0005069



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against human TNF- α (cachectin, TNFSF2). Designed for intracellular flow cytometry. Recognizes human TNF- α ; ideal for inflammation, autoimmune, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human TNF- α mAb
Remarks/Alias	DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha
Species	Human
GeneID (Human)	7124
GeneID	7124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-233 of human TNF- α (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01375
Protein Weight	26kDa
Shipping	Ice bag