

PE/Cyanine7 Rabbit Anti-Human/Monkey Cd19 Mab - Cd19

Item Number: CM0004649

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting CD19 for flow cytometry. Cross-reacts with human and cynomolgus monkey. Recombinant immunogen, 63kDa target. Ideal for B-cell lymphoma and immunotherapy research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human/Monkey CD19 mAb
Remarks/Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID (Monkey)	102145514
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19 (XP_005591597.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Protein Weight	63kDa
Shipping	Ice bag