


PE/Cyanine7 Rabbit Anti-Human/Monkey Cd20 Monoclonal Antibody For Fc - P11836

Item Number: CM0005912



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD20, validated for flow cytometry. Cross-reacts with cynomolgus monkey. Ideal for B-cell development, differentiation, and immunotherapy research applications in immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human/Monkey CD20 mAb
Remarks/Alias	B1; S7; Bp35; CD20; FMC7; CVID5; LEU-16
Species	Human
Gene ID (Human)	931
Gene ID	931
Immunogen	Synthetic Peptide corresponding to a sequence within amino acids 100-200 of human CD20 mAb (NP_068769.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P11836
Protein Weight	14kDa/33kDa
Shipping	Ice bag