

PE/Cyanine7 Rabbit Anti-Human/Monkey Cd25 Monoclonal Antibody For Flow Cytometry - P01589

Item Number: CM0008215



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human and cynomolgus monkey CD25. Designed for flow cytometry detection of IL-2 receptor alpha on regulatory and activated T cells. Research use only.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human/Monkey CD25 mAb
Remarks/Alias	p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Human
Gene ID (Human)	3559
Gene ID	3559
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-213 of human CD25 (NP_000408.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01589
Protein Weight	31kDa
Shipping	Ice bag