

Pe/Cyanine7 Rabbit Anti-Human/Monkey Cd3 Monoclonal Antibody For Fc - P07766

Item Number: CM0008393



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD3 epsilon chain. Validated for flow cytometry applications. Ideal for T-lymphocyte identification and immune research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human/Monkey CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
Gene ID (Human)	916
Gene ID	916
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human CD3 (NP_000724.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag