

Pe/Cyanine7 Rabbit Anti-Human/Monkey Cd45Ra Mab For Fc - P08575-8

Item Number: CM0010571

Introduction

PE/Cyanine7 rabbit monoclonal antibody against CD45RA (P08575-8 isoform). Validated for human and cynomolgus cross-reactivity, suitable for flow cytometry detection of the 138 kDa leukocyte antigen.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human/Monkey CD45RA mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-482 of human CD45RA (P08575-8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag