


PE/Cyanine7 Rabbit Anti-Human/Monkey Cd8A Monoclonal Antibody For Flow Cytometry - P01732

Item Number: CM0008197



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against human and cynomolgus monkey CD8a (SWISS: P01732). Validated for flow cytometry (FC) to detect CD8+ T-cell and NK cell subsets.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human/Monkey CD8a mAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha
Species	Human
GeneID (Human)	925
GeneID	925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8A (NP_001759.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01732
Protein Weight	21kDa/25kDa/30kDa
Shipping	Ice bag