

## PE/Cyanine7 Rabbit Anti-Human/Monkey IgD Monoclonal Antibody For Fc - P01880

Item Number: CM0011960



### Introduction

PE/Cyanine7 Rabbit anti-Human/Monkey IgD Monoclonal Antibody for flow cytometry. Targets IgD (Immunoglobulin heavy constant delta) for B-cell studies and immune signaling research. Human and cynomolgus reactivity. P01880.

[Learn More](#)

| Parameter        | Value                                                                                             |
|------------------|---------------------------------------------------------------------------------------------------|
| Product Name     | PE/Cyanine7 Rabbit anti-Human/Monkey IgD mAb                                                      |
| Remarks/Alias    | Immunoglobulin heavy constant delta                                                               |
| Species          | Human                                                                                             |
| Gene ID          | 3495                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD |
| Source           | Rabbit                                                                                            |
| Category         | Monoclonal Antibodies                                                                             |
| Application      | FC                                                                                                |
| Cross Reactivity | Human, Cynomolgus                                                                                 |
| SWISS            | P01880                                                                                            |
| Protein Weight   | 47kDa                                                                                             |
| Shipping         | Ice bag                                                                                           |