

PE/Cyanine7 Rabbit Anti-Human/Monkey Pd-1/Cd279 Monoclonal Antibody For Fc - Q15116

Item Number: CM0012785



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human and monkey PD-1/CD279. Validated for flow cytometry, this immune checkpoint antibody is ideal for studying T-cell exhaustion and cancer immunotherapy applications.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human/Monkey PD-1/CD279 mAb
Remarks/Alias	PD1; PD-1; CD279; SLEB2; hPD-1; hPD-L; hSLE1
Species	Human
GeneID (Human)	5133
GeneID	5133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-168 of human PD-1/CD279 (NP_005009.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q15116
Protein Weight	32kDa
Shipping	Ice bag