

PE/Cyanine7 Rabbit Anti-Mouse Cd105/Endoglin Monoclonal Antibody For Flow Cytometry - Q63961

Item Number: CM0008596

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD105/Endoglin. Validated for flow cytometry, it enables detection of endothelial cells and is suitable for angiogenesis, cardiovascular, and TGF-beta signaling research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD105/Endoglin mAb
Remarks/Alias	Endo; CD105; S-endoglin
Species	Mouse
Gene ID (Human)	2022
Gene ID (Mouse)	13805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-581 of mouse CD105 (NP_031958.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q63961
Protein Weight	70kDa
Shipping	Ice bag