

PE/Cyanine7 Rabbit Anti-Mouse Cd150/Slam Monoclonal Antibody For Flow Cytometry - Q9Qum4

Item Number: CM0008376

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD150 (SLAM). For flow cytometry, it detects the self-ligand receptor modulating immune cell activation and differentiation. Suitable for immunophenotyping and research on SLAM-mediated signaling in T cells, B cells, and macrophages.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD150/SLAM mAb
Remarks/Alias	Slam, CD150, IPO-3, CDw150, ESTM51, 4933415F16
Species	Mouse
Gene ID (Human)	6504
Gene ID	27218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-242 of mouse CD150/SLAM (NP_038758.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QUM4
Protein Weight	36kDa/38kDa
Shipping	Ice bag