

Pe/Cyanine7 Rabbit Anti-Mouse Cd19 Monoclonal Antibody For Flow Cytometry - P25918

Item Number: CM0010373

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD19, validated for flow cytometry. Targets B-lymphocyte antigen CD19 (UniProt P25918), a B-cell coreceptor involved in BCR signaling and immunity.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD19 mAb
Species	Mouse
GeneID (Human)	930
GeneID	12478
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 19-287 of mouse CD19 (NP_033974.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag