


PE/Cyanine7 Rabbit Anti-Mouse Cd206/Mmr Monoclonal Antibody For Flow Cytometry - Q61830

Item Number: CM0012285



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD206/MMR, ideal for flow cytometry detection of macrophages. Validated on mouse samples. Part of CamelBio's one-stop IVD raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD206/MMR mAb
Remarks/Alias	MR; CD206
Species	Mouse
GeneID (Human)	4360
GeneID	17533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-154 of mouse CD206/MMR (NP_032651.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra), FC
Cross Reactivity	Mouse
SWISS	Q61830
Protein Weight	165 kDa
Shipping	Ice bag