

PE/Cyanine7 Rabbit Anti-Mouse Cd223/Lag-3 Monoclonal Antibody For Flow Cytometry - Q61790

Item Number: CM0004877



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against mouse LAG-3 (CD223), validated for flow cytometry. Suitable for immune checkpoint regulation, T-cell exhaustion, and tumor immunology studies in mouse models.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD223/LAG-3 mAb
Remarks/Alias	Ly66; CD223; LAG-3
Species	Mouse
GeneID (Human)	3902
GeneID	16768
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-442 of mouse CD223/LAG-3 (NP_032505.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61790
Protein Weight	57 kDa
Shipping	Ice bag