

PE/Cyanine7 Rabbit Anti-Mouse Cd23 Monoclonal Antibody For Flow Cytometry - P20693

Item Number: CM0009490

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD23 (FcεRII), validated for flow cytometry. Low-affinity IgE receptor essential for B cell differentiation and IgE regulation in murine immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD23 mAb
Remarks/Alias	CD23; Fce2; Fcer2; Ly-42
Species	Mouse
GenelD (Human)	2208
GenelD	14128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-331 of mouse CD23 (NP_038545.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P20693
Protein Weight	38kDa
Shipping	Ice bag