

Pe/Cyanine7 Rabbit Anti-Mouse Cd25 Monoclonal Antibody For Fc - P01590

Item Number: CM0012637

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD25 (IL-2RA). Optimized for flow cytometry applications. Recognizes the 31 kDa interleukin-2 receptor alpha subunit, key in Treg-mediated immune tolerance.

[Learn More](#)

Parameter	Specification
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43
Species	Mouse
GeneID (Human)	3559
GeneID	16184
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01590
Protein Weight	31kDa
Shipping	Ice bag