

PE/Cyanine7 Rabbit Anti-Mouse Cd34 Monoclonal Antibody For Fc - Q64314

Item Number: CM0009139

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD34 for flow cytometry. Reacts with mouse. Suitable for hematopoietic stem cell identification, bone marrow analysis, and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Mouse
Gene ID (Human)	947
Gene ID	12490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-287 of mouse CD34 (NP_598415.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64314
Protein Weight	41kDa
Shipping	Ice bag