

PE/Cyanine7 Rabbit Anti-Mouse Cd38 Monoclonal Antibody For Flow Cytometry - P56528

Item Number: CM0011678



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD38 for flow cytometry. Reacts with mouse CD38, a transmembrane glycoprotein involved in calcium mobilization and immune cell signaling. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD38 mAb
Remarks/Alias	I-19; ADPRC 1; Cd38-rs1
Species	Mouse
Gene ID (Human)	952
Gene ID	12494
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P56528
Protein Weight	34kDa
Shipping	Ice bag