

Pe/Cyanine7 Rabbit Anti-Mouse Cd4 Monoclonal Antibody For Flow Cytometry - P06332

Item Number: CM0006148

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD4. Validated for flow cytometry. Ideal for studying T-helper cell development, immune response, and macrophage differentiation.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
Gene ID (Human)	920
Gene ID	12504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag