

PE/Cyanine7 Rabbit Anti-Mouse Cd48 Monoclonal Antibody (Flow Cytometry) - P18181

Item Number: CM0008031



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD48 (SLAMF2), validated for flow cytometry. Ideal for studying lymphocyte activation, NK cell function, and immune synapse formation.

[Learn More](#)

Parameter	Specification
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD48 mAb
Remarks/Alias	BCM1; BLAST; Bcm-1; BLAST1; SLAMF2; Sgp-60; BLAST-1; MEM-102
Species	Mouse
Gene ID (Human)	962
Gene ID	12506
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-216 of mouse CD48 (NP_031675.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18181
Protein Weight	27kDa
Shipping	Ice bag