

Pe/Cyanine7 Rabbit Anti-Mouse Cd56/Ncam Monoclonal Antibody For Flow Cytometry - P13595

Item Number: CM0007778



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD56/NCAM. Validated for flow cytometry, enabling precise detection of neural cell adhesion molecule in mouse samples. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD56/NCAM mAb
Remarks/Alias	CD56; Ncam; E-NCAM; NCAM-1
Species	Human
Gene ID (Human)	4684
Gene ID	17967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-711 of mouse CD56 (NP_001106675.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13595
Protein Weight	119 kDa
Shipping	Ice bag