

PE/Cyanine7 Rabbit Anti-Mouse CD62L Monoclonal Antibody For Flow Cytometry - P18337

Item Number: CM0012439

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD62L (L-selectin) validated for flow cytometry. Binds a key leukocyte adhesion molecule mediating lymphocyte homing to lymph nodes. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD62L mAb
Remarks/Alias	Lnhr; CD62L; LAM-1; Ly-22; Lyam1; LECAM1; Ly-m22; Lyam-1; LECAM-1; L-selectin
Species	Mouse
GeneID (human)	6402
GeneID (mouse)	20343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of mouse CD62L (NP_035476.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18337
Protein Weight	42kDa
Shipping	Ice bag