


PE/Cyanine7 Rabbit Anti-Mouse Cd8A Monoclonal Antibody For Fc - P01731

Item Number: CM0011857



Introduction

PE/Cyanine7-conjugated Rabbit anti-Mouse CD8a monoclonal antibody validated for flow cytometry. Specific for mouse CD8a, a coreceptor for MHC class I, enabling cytotoxic T-cell analysis in immune research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
GeneID (Human)	941
GeneID	12525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag