

Pe/Cyanine7 Rabbit Anti-Mouse Cd90.2 Mab For Fc - P01831

Item Number: CM0012694

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD90.2 (Thy-1) validated for flow cytometry. Suitable for immunophenotyping, T cell marker studies, and investigation of cell-cell interactions in synaptogenesis.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD90.2 mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
Gene ID (Human)	7070
Gene ID (Mouse)	21838
Immunogen	Recombinant protein of mouse CD90.2
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01831
Protein Weight	18kDa
Shipping	Ice bag