


PE/Cyanine7 Rabbit Anti-Mouse Cd95/Fas Monoclonal Antibody For Fc - P25446

Item Number: CM0009973



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD95/FAS (UniProt P25446), validated for flow cytometry. Targets the apoptosis-mediating receptor implicated in T-cell tolerance. Ideal for studies of cell death and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD95/FAS mAb
Remarks/Alias	lpr; APO1; APT1; CD95; TNFR6; Tnfrsf6
Species	Mouse
GeneID (human)	355
GeneID	14102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-169 of mouse CD95/FAS (NP_032013.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25446
Protein Weight	37KDa
Shipping	Ice bag