

PE/Cyanine7 Rabbit Anti-Mouse Ki67 Monoclonal Antibody For Flow Cytometry (Intracellular) - E9Pvx6

Item Number: CM0012019

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse Ki67 (E9PVX6). Designed for intracellular flow cytometry, this reagent marks proliferating cells in mouse samples, aiding cell cycle and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse Ki67 mAb
Remarks/Alias	Ki67; Ki-67; D630048A14Rik
Species	Mouse
GeneID (Human)	4288
GeneID	17345
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1816-2163 of mouse Ki67 (NP_001074586.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	E9PVX6
Protein Weight	351 kDa
Shipping	Ice bag