

PE/Cyanine7 Rabbit Anti-Mouse/Rat Cd27 Mab - P41272

Item Number: CM0010544

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse/rat CD27 (TNFRSF7) for flow cytometry. Validated for mouse and rat samples. Costimulatory immune-checkpoint receptor on T, NK, and B cells.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse/Rat CD27 mAb
Remarks/Alias	S152; Tp55; Tnfrsf7
Species	Mouse
GeneID (Human)	939
GeneID	21940
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-182 of mouse CD27 (NP_001028298.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	P41272
Protein Weight	28kDa
Shipping	Ice bag