

Pe/Cyanine7 Rabbit Anti-Mouse/Rat Cd31/Pecam1 Monoclonal Antibody For Flow Cytometry - Q08481

Item Number: CM0005026

Introduction

PE/Cyanine7-conjugated rabbit anti-mouse/rat CD31 (PECAM1) monoclonal antibody for flow cytometry. Targets leukocyte transendothelial migration and cell adhesion. Validated for vascular biology and inflammation research, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse/Rat CD31/PECAM1 mAb
Remarks/Alias	Cd31; Pecam; PECAM-1
Species	Mouse
Gene ID (Human)	5175
Gene ID	18613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-220 of mouse CD31 (NP_001292086.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Mouse, Rat
SWISS	Q08481
Protein Weight	81kDa
Shipping	Ice bag