

Pe/Cyanine7 Rabbit Anti-Nhp Igd Mab - Cynomolgus Monkey Igd

Item Number: CM0010215

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against non-human primate IgD, validated for flow cytometry in cynomolgus and rhesus monkey samples. Ideal for immunological studies and diagnostic assay development.

[Learn More](#)

Parameter	Details
Product Name	PE/Cyanine7 Rabbit anti-NHP IgD mAb
Species	Cynomolgus monkey
Immunogen	Recombinant protein of cynomolgus monkey immunoglobulin delta heavy chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Cynomolgus, Rhesus
Shipping	Ice bag