

Pe/Cyanine7 Rat Anti-Mouse Cd3 Monoclonal Antibody For Flow Cytometry - P04235 / P22646

Item Number: CM0005399

Introduction

PE/Cyanine7-conjugated rat monoclonal antibody targeting mouse CD3 for flow cytometry. Recognizes CD3 delta and epsilon chains, validated for mouse samples. Ideal for T-cell studies.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rat anti-Mouse CD3 mAb
Remarks/Alias	T3d; CD3; T3e; CD3epsilon
Species	Mouse
GenelD (Human)	916/915
GenelD	12500/12501
Immunogen	Synthetic peptide derived from mouse CD3
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P04235/P22646
Protein Weight	19kDa/21kDa
Shipping	Ice bag