

Per2 Rabbit Monoclonal Antibody (Human) For Wb, Elisa, Chip - O15055

Item Number: CM0004078



Introduction

PER2 rabbit recombinant monoclonal antibody validated for WB, ELISA, and ChIP. Specific to human PER2 (O15055), a core circadian clock regulator. Ideal for circadian rhythm and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	PER2 Rabbit mAb
Remarks/Alias	FASPS; FASPS1; PER2
Species	Human
GeneID (Human)	8864
GeneID	8864
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PER2 (O15055).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human
SWISS	O15055
Protein Weight	137 kDa
Shipping	Ice bag