

Percp Anti-Cd34 Monoclonal Antibody For Flow Cytometry - P28906

Item Number: CM0009888

Introduction

PerCP-conjugated rabbit monoclonal antibody targeting human CD34, validated for flow cytometry. Recognizes hematopoietic progenitor cell antigen; suitable for stem cell and endothelial research.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Human CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Human
GeneID (Human)	947
GeneID	947
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-290 of human CD34 (NP_001020280.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28906
Protein Weight	41kDa
Shipping	Ice bag