

Percp Anti-Foxp3 Monoclonal Antibody For Flow Cytometry (Intracellular) - Q9Bzs1

Item Number: CM0007769



Introduction

PerCP-conjugated rabbit monoclonal antibody targeting human/mouse/monkey FOXP3, validated for intracellular flow cytometry. Ideal for regulatory T-cell research, immune tolerance studies. Reacts with human, mouse, cynomolgus. ~47 kDa.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Human/Mouse/Monkey FOXP3 mAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER
Species	Human
GeneID	50943
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXP3 (NP_054728.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Mouse, Cynomolgus
SWISS	Q9BZS1
Protein Weight	47 kDa
Shipping	Ice bag