

Percp Anti-Mouse Cd45 (Ptprc) Monoclonal Antibody For Flow Cytometry - P06800

Item Number: CM0004486

Introduction

PerCP-conjugated rabbit monoclonal antibody targeting mouse CD45 (Ptprc) for flow cytometry (FC). Validated in mouse, detects 130-144 kDa isoforms. Ideal for immune cell profiling.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Mouse CD45 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GeneID (Human)	5788
GeneID	19264
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag