

Percp-Conjugated Rabbit Anti-Mouse Cd4 Monoclonal Antibody For Fc - P06332

Item Number: CM0010056

Introduction

PerCP-conjugated rabbit monoclonal antibody against mouse CD4, validated for flow cytometry (FC). Targets the T-cell surface glycoprotein CD4 (SWISS P06332), ideal for immune response studies and T-helper cell analysis.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
Gene ID (Human)	920
Gene ID	12504
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag