

Percp Mouse Anti-Human Cd3 Monoclonal Antibody For Fc - P07766

Item Number: CM0012370



Introduction

PerCP Mouse anti-Human CD3 Monoclonal Antibody, for flow cytometry, targeting the CD3 epsilon chain on T cells. Ideal for T-cell immunophenotyping and adaptive immunity research. - P07766

[Learn More](#)

Parameter	Value
Product Name	PerCP Mouse anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GeneID (Human)	916
GeneID	916
Immunogen	Recombinant protein of human CD3
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag