

Percp Mouse Igg 2A Isotype Control - Isotype Control

Item Number: CM0010775

Introduction

High-quality PerCP-conjugated mouse IgG 2a isotype control monoclonal antibody for flow cytometry. Ideal for non-specific binding assessment and accurate gating in FC assays. Reliable IVD raw material from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	PerCP Mouse IgG 2a isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag