

Percp Mouse Igg 2A Isotype Control

Item Number: CM0012532

Introduction

PerCP-conjugated Mouse IgG2a isotype control monoclonal antibody for flow cytometry. Sourced from mouse, this reagent serves as a negative control to evaluate non-specific binding in FC assays, ensuring accurate data interpretation.

[Learn More](#)

Parameter	Specification
Product Name	PerCP Mouse IgG 2a isotype control
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag