

Percp Rabbit Anti-Human/Monkey Cd45 Monoclonal Antibody For Flow Cytometry - P08575

Item Number: CM0008816



Introduction

PerCP-conjugated rabbit monoclonal antibody targeting human CD45 (P08575). Validated for flow cytometry in human and cynomolgus samples. Ideal for immunophenotyping, T-cell research, and immune profiling.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575
Protein Weight	131kda/136kDa/141kDa/143kDa/147kDa
Shipping	Ice bag