


Percp Rabbit Anti-Human/Monkey Cd8A Mab - P01732

Item Number: CM0008997



Introduction

PerCP-conjugated rabbit anti-human/monkey CD8a monoclonal antibody for flow cytometry (FC). Recognizes CD8 alpha chain (P01732), validated for human and cynomolgus monkey samples. Ideal for immunological research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Human/Monkey CD8a mAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha
Species	Human
GeneID (Human)	925
GeneID	925
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8A (NP_001759.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01732
Protein Weight	21kDa/25kDa/30kDa
Shipping	Ice bag