


Percp Rabbit Anti-Mouse Cd209B/Dc-Signr1 Monoclonal Antibody For Fc - Q8Cj91

Item Number: CM0008145



Introduction

PerCP-conjugated rabbit monoclonal antibody against mouse CD209b/DC-SIGNR1 (SIGNR1), validated for flow cytometry. Reacts with mouse. Ideal for dendritic cell and pathogen receptor studies.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Mouse CD209b/DC-SIGNR1 mAb
Remarks/Alias	OtB7; SIGNR1; mSIGNR1; DC-SIGNR1; 1810030I22Rik
Species	Mouse
GeneID (Human)	30835
GeneID	69165
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8CJ91
Protein Weight	16kDa/33kDa/36kDa/37kDa
Shipping	Ice bag