


Percp Rabbit Anti-Mouse Cd44 Monoclonal Antibody For Fc - P15379

Item Number: CM0005760



Introduction

PerCP-conjugated rabbit monoclonal antibody against mouse CD44 (Ly-24, Pgp-1, HERMES). Validated for flow cytometry, suitable for lymphocyte homing, adhesion, and activation studies.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Mouse CD44 mAb
Remarks/Alias	Ly-24, Pgp-1, HERMES
Species	Mouse
Gene ID (Mouse)	12505
Gene ID (Human)	960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	86 kDa/81 kDa/76 kDa/72 kDa
Shipping	Ice bag