


Percp Rabbit Anti-Mouse Cd8A Monoclonal Antibody For Flow Cytometry - P01731

Item Number: CM0011714



Introduction

PerCP-conjugated rabbit monoclonal antibody against mouse CD8a (Lyt-2) for flow cytometry. Validated for mouse samples, recognizes the T-cell coreceptor involved in MHC class I interaction and CTL activation.

[Learn More](#)

Parameter	Value
Product Name	PerCP Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
GeneID (human)	925
GeneID	12525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag