

Percp Rabbit Igg Isotype Control

Item Number: CM0010963

Introduction

PerCP-conjugated Rabbit IgG isotype control monoclonal antibody designed for flow cytometry (FC) applications. This rabbit-derived antibody serves as an essential negative control to assess nonspecific binding, ensuring accurate gating and data interpretation in immunological assays.

[Learn More](#)

Parameter	Specification
Product Name	PerCP Rabbit IgG isotype control
Immunogen	A synthesized peptide derived from rabbit IgG isotype control
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag