

Percp/Cyanine5.5 Anti-Mouse Cd11C Monoclonal Antibody For Fc - Q9Qxh4

Item Number: CM0011078

Introduction

PerCP/Cyanine5.5-conjugated anti-mouse CD11c monoclonal antibody (clone N418, Armenian hamster) for flow cytometry. Validated in mouse. Research use for inflammatory response and monocyte adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 anti-Mouse CD11c mAb
Remarks/Alias	Cr4; N418; Cd11c
Species	Mouse
GeneID (human)	3687
GeneID	16411
Immunogen	Recombinant protein of mouse CD11c
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QXH4
Protein Weight	129kDa
Shipping	Ice bag